

Competitive Product Test Results (Surat Thani 7/02)

GROUP	BRAND NAME	BATCH NO.	DIMENSIONS						BURST PROPERTIES				LUBRICANT QUANTITY		SEAL INTEGRITY	
			Length (mm)		Width (mm)		Thickness (mm)		Sample size	MV-1	NCV-2	MP-3	NCP-4	Avg.	Range	NC
			Avg.	Range	Avg.	Range	Avg.	Range								
		SPEC.							Min. 18.0		Min. 1.0					
Lubricated	Trojan-enz Lubricated	093750	181	178-185	53.5	53.5-54.0	0.076	0.071-0.080	200	41.94	0	1.96	0	701	627-739	0
	Durex Extra Sensitive	013000	181	180-185	52.5	52.0-52.5	0.058	0.053-0.063	200	35.60	0	1.56	0	426	416-443	0
	Okamoto	370012	185	183-188	54.0	53.5-54.5	0.055	0.051-0.066	200	35.66	2	1.66	3	140	118-179	0
	LifeStyles (1512)	0206090616	184	181-185	53.5	53.5-54.0	0.070	0.068-0.072	200	41.84	0	1.57	0	476	414-515	0
Non- Lubricated	Trojan Non-Lubricated	090750	187	184-190	53.5	53.0-535	0.071	0.067-0.077	200	35.70	1	2.10	0	Non - Lubricated		3
	Durex Natural Feeling Non-Lub	002800	184	180-187	52.0	51.5-52.0	0.066	0.063-0.070	200	40.48	1	2.02	1	Non - Lubricated		0
	LifeStyles (4412)	0206102416	185	180-189	54.0	53.5-54.0	0.070	0.068-0.071	200	43.82	0	1.71	3	Non - Lubricated		0
Spermicidal	Trojan-enz Spermicidal	093250	184	178-194	53.5	53.0-54.0	0.074	0.069-0.080	200	42.56	0	1.68	0	584	525-619	0
	Durex Extra Sensitive	012800	182	175-187	52.0	51.5-52.5	0.064	0.063-0.070	200	35.20	2	1.64	2	552	521-575	1
	Okamoto	390012	188	185-190	54.0	53.5-54.5	0.054	0.050-0.057	200	34.80	0	1.47	1	541	536-547	0
	LifeStyles (1812)	0111182416	189	185-191	54.0	53.5-54.0	0.071	0.069-0.072	200	38.95	2	1.60	3	353	338-371	0
Thin (Lubricated)	Trojan Ultra Thin	092640	179	174-182	53.5	53.0-54.0	0.071	0.063-0.074	200	34.68	1	1.66	1	366	320-412	4
	Durex Ultimate Feeling	024900	189	182-195	52.0	51.0-53.5	0.064	0.057-0.073	198	31.76	2	1.97	1	512	490-544	5
	LifeStyles (1712)	0203112716	183	180-188	53.0	53.0-54.5	0.068	0.066-0.069	200	41.30	0	1.66	0	508	477-532	0
Non-Latex (Lambskin)	Naturalamb	098750	200	190-216	71.5	69.0-73.5	0.074	0.040-0.135	182	1.48	181	3.26	24	2064	1779-2765	0
Non-Latex (Synthetic)	Trojan Supra	090340	179	170-186	58.0	57.5-59.0	0.055	0.050-0.059	200	14.74	154	2.06	0	758	742-789	0
	Durex Avanti	058000	178	159-183	60.5	60.0-61.0	0.050	0.045-0.055	200	13.69	197	3.72	0	710	656-749	0
X-Large	Trojan Magnum	064212	192	190-195	55.5	55.0-56.0	0.076	0.072-0.082	200	38.14	0	1.66	1	567	524-636	1
	LifeStyles (1112)	0203032416	189	186-190	56.5	56.0-57.0	0.074	0.072-0.075	200	36.60	0	1.62	0	480	403-550	0
Color/Flavors	Durex Colors & Scents	08300	186	184-195	52.0	51.0-52.5	0.069	0.065-0.073	200	38.31	0	2.00	0	440	357-469	1
	LifeStyles (2812)	0108832716	188	181-192	53.0	53.0-54.0	0.070	0.066-0.074	200	36.30	1	1.64	1	491	469-516	0
Desensitizing	Trojan Extended Pleasure	097250	183	180-188	54.0	53.5-54.5	0.081	0.074-0.088	200	44.65	0	1.90	0	863	746-979	0

**-1 = Mean volume score for the lot tested (18.0 is minimum requirement) -2 = Number of condoms which failed to meet 18.0 requirement
-3 = Mean pressure score for the lot tested (1.0 is minimum requirement) -4 = Number of condoms which failed to meet 1.0 requirement**